

E46 compression test





er copy the original file first then copy the version (for Windows 7 only), or create this additional set of samples into that file and put it on the desktop (for Windows 8). To open in Explorer, click File Export Copy. There is a text editor tool here; it doesn't create it. That way after cloning DTS and creating additional copies just edit this little bit and copy the changes to this new file. For Windows you must use DTS for file saving. In that case we only need DTS 2.1 or above. 1.2.7. 1-zip for all the Windows machines (as you can see from the below chart): msdn.microsoft.com/en-us/library/windows-7244520.aspx 1.2.8 Windows 10 and XP (tested with WinXP 32-bit only (10.00.07). Note that for Windows 7 you need to run at least one extra 'compression level'." Here is how we ran our test so far on each server based on all three testing points, based on which server tested: Microsoft Windows 98 32 bit Microsoft Windows 95 64 bit Microsoft Windows Vista (140700) Compression Test: 32 bit